



Multi-Patient Disposable Set

for MEDRAD[®] MRXperion MR Injection System



Enhance patient care



Streamline workflow efficiency



Reduce waste



MR Multi-Patient Disposables	Catalog	Box Quantity
Multi-Patient Kit (24-hour use)	XP MP	10
Single-Patient Line (250cm)	XP SP	50
Clamp (optional)	MRX MP CMP	100

* Fluid colors are for illustrative purposes only.

MR Multi-Patient Disposable Set

Streamlining your workflow with one syringe kit a day for multiple studies.



Tubing Design
250 cm Single Patient Disposable

115 mL syringe

65 mL syringe

Multi-Patient Spike
Spikes can be used to access multiple containers of same contrast and saline for up to 24 hours.

Multi-Guard Connector
A recessed connection prevents accidental contamination of fluid path. The MEDRAD® MRXperion Multi-Patient System includes a 15 cm multi-patient tube that connects to the syringes and is out of patients' reach.

Multi-Patient Disposable Kit Contains:

- > 1 x 65 mL syringe and 1 x 115 mL syringe
- > Multi-Patient Transfer Sets
- > Multi-Patient Connector Tube


Use one kit for up to 24 hours



Single-Patient Line
250cm single patient disposable includes two check valves to prevent backflow from the patient.



Optional Tube Clamp
Clamp will minimize dripping to reduce the risk of contamination.



Bayer reserves the right to modify the specifications and features described herein or to discontinue any product or service identified in this publication at any time without prior notice or obligation. Please contact your authorized Bayer representative for the most current information.

Bayer, the Bayer Cross, MEDRAD, MEDRAD MRXperion, and MRXperion are trademarks owned by and/or registered to Bayer in the U.S. and/or other countries. Other trademarks and company names mentioned herein are properties of their respective owners and are used herein solely for informational purposes. No relationship or endorsement should be inferred or implied. The individuals depicted in this presentation are actors and not actual health care providers or patients.

If you want to report a side effect or quality complaint, please contact your healthcare professional (e.g., physician or pharmacist) or your local Health Authority. Reports can also be directed to <https://safetrack-public.bayer.com/>.

SRV004E © 2025 Bayer. This material may not be reproduced, displayed, modified or distributed without the express prior written consent of Bayer.



Bayer Inc.
2920 Matheson Blvd East
Mississauga (Ontario) L4W 5R6
Telephone: 800-268-1432
Fax: 800-567-1710