



# MEDRAD<sup>®</sup> Mark 7 Arterion Injection System



**MEDRAD<sup>®</sup>** Mark 7 Arterion  
Injection System

MEDRAD® Mark 7 Arterion has stylish ergonomics that are designed for function; the pedestal with a small footprint and reduced weight increases mobility around your busy lab. The head, now lighter, is designed with an integrated handle to make it easier to position for injection.



\*As compared with MEDRAD® Mark V ProVis Angiographic Injection System

\*\* Refer to the Operator Manual for a description of the intended use, contraindications, warnings and instructions of use.

# Designed to Support Your Workflow

## Outstanding Syringe Clarity

The clear polycarbonate barrel of the MEDRAD® Mark 7 Arterion Syringe may enable improved visualization of both contrast and air, facilitating your safety monitoring.

### Color Touchscreen Display

- Injection protocol, contrast media volume in syringe and injection results are displayed on the main screen
- Store and sort up to 35 injection protocols
- Review history of up to 50 actual injections

### Ergonomic Design

- Handle on injector head
- Large handle on pedestal for mobility
- Large, smooth rolling casters lock in place
- 20% smaller footprint and 32% lighter compared to MEDRAD® Mark V ProVis Injection System

### Articulating Front Load Syringe

- Easy set up, use and tear down

### Clinical Application Service

- On-site clinical training to help your staff perform procedures with confidence from day one



# MEDRAD® Mark 7 Arterion Injection System Ordering Information

Product	Catalog Number
<b>MEDRAD® Mark 7 Arterion Injection System</b>	
Pedestal	ART 700 PEDL
Table Mount	ART 700 TABL
Freestanding Pedestal on Wheels*	KMA 320RT
Adjustable Height Pedestal on Wheels*	KMA 330
Dual Display Options*	ART 700 2DCU
Dual Display Adjustable Table Mount	A-DD-ADJ-TM
Dual Display Fixed Table Mount	A-DD-FIX-TM
Dual Display Wall Mount	A-DD-WALL
Dual Display Desk Mount	A-DD-DESK
Foot Switch Option*	ART 700 FS
<b>Sterile Disposables</b>	
Disposable Syringe with Quickfill Tube, 150 mL (25/box)*	ART 700 SYR

\* Optional



2920 Matheson Blvd East  
Mississauga ON L4W 5R6  
Phone: (800) 268-1432  
Fax: (800) 567-1710

Bayer reserves the right to modify the specifications and features described herein or to discontinue any product or service identified in this publication at any time without prior notice or obligation. Please contact your authorised Bayer representative for the most current information.

All patient data appear in this document are fictitious. No actual patient information is shown.

Bayer Cross, MEDRAD, P3T, Virtual Care and Radimetrics are trademarks owned by and/or registered to Bayer in the U.S. and /or other countries. Other trademarks and company names mentioned herein are properties of their respective owners and are used herein solely for informational purposes. No relationship or endorsement should be inferred or implied.

©2023 Bayer. This material may not be reproduced, displayed, modified or distributed without the express prior written consent of Bayer.