



The next big step in CT injector technology

CEntargo CT Injection System

What is MEDRAD[®] Centargo?

The latest addition to Bayer's range of CT products is MEDRAD[®] Centargo, an innovation CT injection system that simplifies workflow and provides efficiencies for radiology departments.

The workload of radiology departments is continuously increasing. This growing demand is putting pressure on staff and equipment, leading to an increasing need to do more with less. Centargo has been designed to maximize efficiency and automate workflow, allowing radiology staff to spend more time with their patients.

Efficient daily set-up The daily set-up process can be completed in under 2 minutes.

Loading Day Set

Loading Connecting Contrast Patient

line

Quick patient-to-patient set-up The simple snap-in patient line auto-primes upon insertion and is ready for the next patient in less than 20 seconds.

and Saline



Dual screen configuration

Touch screen displays both in the control room and the scan room so that staff can be close to their patient.



Automated documentation Integrated barcode reader reduces manual data entry and

provides easy traceability and access to important contrast and injection details



Flexible configurations Cordless mobility through battery power and a Wi-Fi enabled configuration.



The value of MEDRAD[®] Centargo for radiology staff and patients



Confidence Efficiency Quality Centargo is designed for a simple, Less manual data entry Centargo offers trusted Bayer reliable and automated workflow, easy and an intuitive user MEDRAD[®] quality, relying on decades of experience, daily set-up and quick interface help staff patient-to-patient changeover confidently complete comprehensive service between scans: their daily daily tasks: networks and expert clinical support Daily set-up process that can Touch screen within the scan Inlet and outlet air be completed in under room so that staff can be close detection, in addition to 2 minutes to their patient automatic air removal within reservoirs, to help Integrated barcode reader that prevent air injection **Pre-assembled** reduces manual data entry and 24-hour Day Set provides easy traceability and Piston-based fluid delivery and 胍 access to contrast and injection heat maintainers for contrast details Simple snap-in Patient reservoirs de-signed to provide Line that auto-primes quality fluid delivery, even Two check vavles in the upon insertion, and is during the most demanding patient line to prevent ready for the next patient protocols cross-contamination in less than 20 seconds Informatics features that >-Advanced scanner connectivity, automate documentation and R synchronizing scan timing and calculate individualized simplifying workflow* injection protocols Cordless mobility with battery power and a Wi-Fi enabled configuration

* requires additional purchase of ISI2



MEDRAD® Centargo - Allowing radiology staff to do less and care more



© 2022 Copyright Bayer AG. Bayer, the Bayer Cross, MEDRAD, MEDRAD Centargo, Centargo, P3T and Radimetrics are trademarks of the Bayer group of companies. All rights reserved.